

Monitoring cell membrane depolarization due to L12 electroporation using fluorescent plasma membrane potential indicator

Anja Blažič, Lea Rems

University of Ljubljana, Faculty of Electrical Engineering

Duration of the experiments: 90 min

Max. number of participants: 4

Location: Cell Culture Laboratory 3

Level: Advanced

PREREQUISITES

Participants should be familiar with Laboratory safety (S1). The basic knowledge of handling with cells is required for this laboratory practice.

BACKGROUND

All cells maintain an electric potential difference across their plasma membranes, which results from the differences in membrane permeabilities for potassium, sodium, calcium, and chloride ions. This potential difference is called the resting transmembrane voltage or resting potential and is maintained by a system of ion channels and pumps. By convention the resting transmembrane voltage is negative, meaning that the cell interior is electrically more negative compared to its exterior. The value of the resting transmembrane voltage changes dynamically with the cell cycle and has an important biological function by controlling the activity of various membrane proteins. When cells are electroporated, their transmembrane voltage changes and the membrane remains depolarized for several minutes after pulse exposure. As membrane depolarization acts as a biological signal, factors that influence prolonged depolarization upon electroporation can have an important influence on the biological outcome of electroporation.

The aim of this laboratory practice is to monitor the time course of transmembrane voltage changes after exposure to conventional electroporation pulses, and to determine the influence of pulse amplitude and ambient temperature on the extent and longevity of membrane depolarization.

EXPERIMENT

We will monitor the time course of membrane potential changes using the plasma membrane potential indicator (PMPI) of the FLIPR Membrane Potential Assays Kit (Molecular Devices). PMPI consists of a two-part system which includes a fluorescent anionic voltage-sensor and a quencher. When the interior of the cell has a relatively negative charge the anion dye remains on the extracellular side, where the quencher prevents fluorescence excitation. During depolarization the voltage sensor translocates to the intracellular side, which increases the fluorescence intensity inside the cells. This translocation is reversible, which makes it possible to monitor membrane depolarization and repolarization as the cells recover and return to their baseline level. The fluorescence intensity is linearly proportional to the change in the membrane potential difference, making this method comparable to patch-clamp measurements upon dye calibration.

Protocol: U87-MG human glioblastoma cells will be plated in Nunc Lab-Tek II chambered coverglass three to four days before the experiments. On the day of the experiments, the cells will be stained for 30 min at 37°C with the Component A of the FLIPR Membrane Potential Assay Red (Molecular Devices, #R7291), diluted in Live Cell Imaging Solution (Invitrogen, #A14291DJ). Afterwards, the cells will be placed on the microscope stage. The cells will be exposed to a single 100 μ s pulse of selected amplitude, delivered by β tech B10 pulse generator through a pair of Pt-Ir wire electrodes. Time lapse images of the cells will be acquired before and after pulse application. We will compare the response of the cells at room temperature and at 37°C, to observe how cell recovery depends on the temperature. The cells will be imaged on inverted microscope Leica DMI8 with LED8 illumination source controlled by the LasX software (all Leica Microsystems). The membrane potential dye will be excited with green LED (555/28 nm) and its fluorescence will be passed through DFT51010 filter and detected with the Leica DFC9000 Gt camera.




Figure 1: Brightfield image (left) and fluorescence (right) of U87-MG cells stained with FLIPR membrane potential assay.

FURTHER READING:

- Wegner L.H., Frey W., Silve A. Electroporation of DC-3F cells is a dual process. *Biophys. J.* 108:1660-71, 2015.
- Baxter D.F., Kirk M., Garcia A.F., Raimondi A., Holmqvist M.H., Flint K.K., Bojanic D., Distefano P.S., Curtis R., Xie Y. A novel membrane potential-sensitive fluorescent dye improves cell-based assays for ion channels. *J. Biomol. Screening* 7:79-85, 2002.
- Klapperstück T., Glanz D., Hanitsch S., Klapperstück M., Markwardt F., Wohlrab J. Calibration procedures for the quantitative determination of membrane potential in human cells using anionic dyes. *Cytometry Part A* 83A:612-26, 2013.
- Burke R.C., Bardet S.M., Carr L., Romanenko S., Arnaud-Cormos D., Leveque P., O'Connor R.P. Nanosecond pulsed electric fields depolarize transmembrane potential via voltage-gated K^+ , Ca^{2+} and TRPM8 channels in U87 glioblastoma cells. *Biochim. Biophys. Acta (BBA) - Biomembranes* 1859:2040-50, 2017.
- Dermol-Černe J., Miklavčič D., Reberšek M., Mekuč P., Bardet S.M., Burke R., Arnaud-Cormos D., Leveque P., O'Connor R.P. Plasma membrane depolarization and permeabilization due to electric pulses in cell lines of different excitability. *Bioelectrochemistry* 122:103-14, 2018.
- Blackiston D.J., McLaughlin K.A., Levin M. Bioelectric controls of cell proliferation: ion channels, membrane voltage and the cell cycle. *Cell Cycle* 8(21):3527-36, 2009.
- Vodovnik L., Miklavčič D., Serša G. Modified cell proliferation due to electrical currents. *Med. Biol. Eng. Comput.* 30:CE21-8, 1992.

NOTES & RESULTS

